Sample preparation for plant tissue

Freeze plant tissues with liquid N₂ Grind plant tissue into powder with a mortar and pestle in liquid N₂. Weighing 300~500 mg (fresh weight) plant tissue and keep the samples in liquid N₂. Add 1.125 mL of a mixture cold CHCl₃/MeOH (2/1) Mix sample for 15 min at 4°C Add 375 µl cold H2O Mix sample for 15 min at 4°C Centrifuge at 13,000 rpm for 10 min at 4°C Move the upper (hydrophilic) layer to new eppendorf Move the lower (hydrophobic) layer to another eppendorf (Pierce the protein disc carefully) Freeze by liquid N₂ Dry out by SpeedVac and keep at -80 °C (before LC-MS analysis) Reconstitute the hydrophilic layer with 180~300 μl MeOH:H2O (1/1) and hydrophobic layer with CHCl3:MeOH (2/1) Centrifuge at 13,000 rpm for 15 min at 4°C and take the supernatant