


User manual for *FeatureHunter* ver.1.4



Release date: Feb, 2024



About *FeatureHunter*



FeatureHunter was developed as the ultimate solution for untargeted nucleic acid adductomics, empowering researchers to effortlessly process data from a wide range of cutting-edge high-resolution mass spectrometry and leading MS vendors (e.g., Agilent, Thermo, and Waters).



Table of the contents

1. System requirements
2. Installation
3. Steps for the data processing within *FeatureHunter*
4. Appendix
 - A1. Installation of Visual C++ redistributable package
 - A2. Installation of .NET Framework
 - A3. Parameter settings for data format conversion in MSconvert
 - A4. Management of Feature List
 - A5. Management of Artificial Adducts List
 - A6. Examples for Tag Usage
 - A7. Configuration of data export from MZmine
 - A8. Scheme of results integration with MZmine




1. System Requirements



Recommended requirements

- 4 or more CPU cores
- 16 GB or more RAM
- 64-bit Windows 10 or later (macOS is not available)



2. Installation

Prerequisites

- Visual C++ redistributable package (64-bit)
- .NET Framework 4.7.2 or later (64-bit)
- FeatureHunter (ver. 1.4)
- MSConvert



- The successful running of *FeatureHunter* requires the installation of Visual C++ redistributable package and .NET Framework.
- A detailed guide for the installation of Visual C++ redistributable package and .NET Framework can be found in [Appendix A1 and A2](#).
- The raw data obtained from vendor format must be converted into .mzML format (please see [Appendix A3](#)) for successfully running in *FeatureHunter*.

Download Software, Feature List and Example Data

Downloads

1

Software

FeatureHunter

2

Feature List

List of 132 tags with descriptions

download

List of tag combinations for ready-to-use

download

3

Example Data

For testing the general function within FeatureHunter

Mixture of modified nucleic acid standards

raw mzML

Mixture of DNA and RNA co-treated with d0- and d8-chlorambucil

raw mzML

For testing the data integration between FeatureHunter and MZmine

Mixture of modified nucleic acid standards

tsv (FeatureHunter) csv (MZmine)

Mixture of DNA and RNA co-treated with d0- and d8-chlorambucil

tsv (FeatureHunter) csv (MZmine)

Aligned QC samples (QC1-QC3)

tsv (FeatureHunter) csv (MZmine)

FeatureHunter can be downloaded from the following link:

https://msomics.abrc.sinica.edu.tw/FeatureHunter/?page_id=9

- 1 Download the **Software**
- 2 Download the **Feature List**, which guides users to choose the tags (features) for extracting compounds of interest
- 3 Download the **Example Data**, available in various formats (Thermo.raw and mzML) and outputs (.tsv and .csv), to freely test and verify the software

Fig-1. The page for downloading the software, the list of tag descriptions and applications, and example data

Checking the download file




Fig-2. The items included in the *FeatureHunter* download folder

Once the download is complete, the folder will contain the following four items:

- ① **bin:** This folder contains the executable files of *FeatureHunter*.
- ② **params:** This folder contains two important parameter files used in *FeatureHunter*: one for feature extraction settings and another for artificial adducts match configurations.
- ③ **FeatureHunter.exe:** The main program executable for *FeatureHunter*.
- ④ **UserManual.pdf:** A tutorial providing a guidance on how to use *FeatureHunter*.



- For successful execution of *FeatureHunter*, place the (1) “bin” folder, (2) “params” folder, and (3) “FeatureHunter.exe” file in the same directory.
- Additionally, it is recommended to store the downloaded folder on the C drive at the path: **C:\FeatureHunter**

3. Steps for the data processing within *FeatureHunter*




Fig-3. The flowchart of the *FeatureHunter* processing

3-1. Launch the FeatureHunter [Step-1]




Fig-4. Starting the software

3-2. Feature extraction (Function-1) [Step-2]

FeatureHunter ver. 1.3

Feature Extraction | Feature Annotation and Classification | Peak Confirmation via MSmine (Optional)

mzML Import:

Browse Add Remove

C:\Mixture_of_DNA_and_RNA_cotreated_with_dO-and_d8-chlorambucil.mzML
C:\Mixture_of_modified_nucleic_acid_standards.mzML

Parameters

Project Name: TEST1

Mass Tolerance: 5 ppm

Ion Mode: ☒ Positive

Feature List: Browse C:\FeatureHunter_v1.3_Official\params\Feature_List.csv

Artificial Adducts: Browse C:\FeatureHunter_v1.3_Official\params\Artificial_Adduct_List.csv

Output Folder: Browse C:\TEST

RUN

① Select the “Feature Extraction”.

② Click “Browse” to import the raw data file in [mzML](#) format.



As the mzML format was required, carefully refer to the parameter settings for data conversion in MSconvert provided in [Appendix A3](#).

Fig-5. Interface of “Feature Extraction” in *FeatureHunter*

3-2. Feature extraction (Function-1) [Step-2]

FeatureHunter ver. 1.3

Feature Extraction | Feature Annotation and Classification | Peak Confirmation via MSmine (Optional)

mzML Import:

C:\Mixture_of_DNA_and_RNA_co-treated_with_dO-and_d8-chlorambucil.mzML
C:\Mixture_of_modified_nucleic_acid_standards.mzML

3 Parameters

Project Name: (a)

Mass Tolerance: ppm (b)

Ion Mode: ☒ Positive (c)

Feature List: C:\FeatureHunter_v1.3_Official\params\Feature_List.csv (d)

Artificial Adducts C:\FeatureHunter_v1.3_Official\params\Artificial_Adduct_List.csv (e)

Output Folder: C:\TEST

Fig-5. Interface of "Feature Extraction" in *FeatureHunter*

- ③ Set the parameters:
- (a) Specify a "Project Name" which will be automatically added as a prefix to the output file name. For example: **TEST1**_file name.
 - (b) Set the "Mass Tolerance" with ppm for the artificial adduct match and the peak height of precursor's EIC search.
 - (c) Positive ion mode is default setting.
 - (d) Import the Feature List file. (The default Feature List.CSV can be found in the params folder)
 - (e) Import the Artificial Adducts List file. (The default Artificial Adducts List.CSV can be found in the params folder)



The Feature List and Artificial Adducts List can be edited and are detailed in **Appendix A4 and A5**, respectively.

3-2. Feature extraction (Function-1) [Step-2]

FeatureHunter ver. 1.3

Feature Extraction | Feature Annotation and Classification | Peak Confirmation via MZmine (Optional)

mzML Import:

Browse Add Remove

C:\Mixture_of_DNA_and_RNA_co-treated_with_dO- and d8-chlorambucil.mzML
C:\Mixture_of_modified_nucleic_acid_standards.mzML

Parameters

Project Name: TEST1

Mass Tolerance: 5 ppm

Ion Mode: ☒ Positive

Feature List: Browse C:\FeatureHunter_v1.3_Official\params\Feature_List.csv

Artificial Adducts Browse C:\FeatureHunter_v1.3_Official\params\Artificial_Adduct_List.csv

4 Output Folder: Browse C:\TEST

5 RUN

④ Specify the output path and folder. (The output files will be in .TSV format.)

⑤ Click "RUN" to start execution.

Fig-5. Interface of "Feature Extraction" in *FeatureHunter*

When processing




Fig-6. Display of the running progress during the execution

Finished




Fig-7. Display of the running completed




Fig-8. The results output in the specified output folder

3-3. Feature Annotation and Classification (Function-2) [Step-3]




Fig-9. Interface of “Feature Annotation and Classification” in *FeatureHunter*

- ① Select the “Feature Annotation and Classification”.
- ② Click “Browse” to import the .TSV file that was generated from Feature Extraction (Function-1).
- ③ Set the parameters:
 - (a) Specify a “Project Name” which will be automatically added as a prefix to the output file name. For example: **TEST2**_TEST1_file name.
 - (b) Set the “Mass Tolerance” with ppm for the feature classification.
 - (c) Set the “RT Tolerance” with seconds for the feature classification.

3-3. Feature Annotation and Classification (Function-2) [Step-3]

FeatureHunter ver. 1.3

Feature Extraction | Feature Annotation and Classification | Peak Confirmation via MZmine (Optional)

TSV Import:

Browse Add Remove

C:\TEST\TEST1_Mixture_of_DNA_and_RNA_co-treated_with_dO-and_d8-chlorambucil.tsv
C:\TEST\TEST1_Mixture_of_modified_nucleic_acid_standards.tsv

Parameters

Project Name: TEST2

Mass Tolerance: 5 ppm

RT Tolerance: 60 seconds

Alignment of the output results

☒ Enable

Mass Tolerance 5 ppm

RT Tolerance: 60 seconds

Tag use in

☒ ∩ (Intersection) 1,2

☐ ∪ (Union)

☐ - (Difference)

Add Clear

Tags Combination [1]∩[2]

Output Folder: Browse C:\TEST

RUN

Fig-9. Interface of “Feature Annotation and Classification” in *FeatureHunter*

- ④ Use tags to classify the annotated features.
- The tags that can be used correspond to the Feature List.csv file located in the params folder.
 - The identity of each tag is described in 'List of 132 tags with descriptions' which can be downloaded from website.
 - Tags can be inputted individually or in combination using the **intersection** (\cap), **union** (\cup), or **difference** ($-$) functions.
 - Multiple tags can be inputted by separating them with a comma.
 - Leaving the field blank is allowed if the function is not being used.



The guide on tag usage can be found in **Appendix A6**.

3-3. Feature Annotation and Classification (Function-2) [Step-3]

FeatureHunter ver. 1.3

Feature Extraction | Feature Annotation and Classification | Peak Confirmation via MZmine (Optional)

TSV Import:

Browse Add Remove

C:\TEST\TEST1_Mixture_of_DNA_and_RNA_co-treated_with_d0-and_d8-chlorambucil.tsv
C:\TEST\TEST1_Mixture_of_modified_nucleic_acid_standards.tsv

Parameters

Project Name: TEST2

Mass Tolerance: 5 ppm

RT Tolerance: 60 seconds

5 Alignment of the output results

☒ Enable

Mass Tolerance: 5 ppm

RT Tolerance: 60 seconds

Tag use in

☐ (Intersection) 1,2

☐ (Union)

☐ (Difference)

Add Clear

Tags Combination [1]∩[2]

Output Folder: Browse C:\TEST

RUN


Fig-9. Interface of “Feature Annotation and Classification” in *FeatureHunter*

⑤ (Optional) Alignment of the output results.

- At least two .TSV file inputs are required for alignment.
- Enabling this option means that the classified results are generated from multiple sample analyses, and they are further aligned to create an aligned feature peak list with precursor ion m/z, retention time (RT), and peak intensity.
- Set the “Mass Tolerance” in ppm to align the ion signals across different datasets.
- Set the “RT Tolerance” in seconds to align the ion signals across different datasets based on their RT.



- Signal alignment was initiated using the m/z and RT values obtained from importing the first dataset (.tsv).
- The notation ‘n/a’ indicated the absence of aligned signals in the current data.



⑥ Specify the output path and folder. (The output files will be in .TSV format.)

⑦ Click "RUN" to start execution.

Fig-9. Interface of “Feature Annotation and Classification” in *FeatureHunter*

Finished




Fig-10. Display of the running completed

Output success




Fig-11. The individual classified peak list (i.e., [1] and [2]) as well as the aligned results, which are output in the specified output folder

An example of output results from feature annotation and classification (Format-1)

■ The output results can be viewed in excel



When using a tag with **category 1 or 2** (see [Appendix A4](#))

Information of detected MS ¹ ions							Information of detected MS ² ions				Information of detected features			
Classification by tag	MS1 Scan Number	RT (min)	Precursor Ion m/z	Precursor Ion Intensity	Precursor Ion Charge	MS Level	Collision Type	Collision Energy	MS2 Scan Number	Product Ion m/z	Product Ion Intensity	Delta Mass (Da)	Mass Error (ppm)	Tags matched
1	2877	3.45	252.108627	1.72E+05	na	2	beam-type collision-induced dissociation	60	2934	136.061523	1.64E+04	116.047104	2.07	1, 124
1	4601	5.43	224.102966	5.57E+05	1	2	beam-type collision-induced dissociation	60	4647	108.055565	1.42E+04	116.047401	0.49	1
1	5897	6.94	294.119659	2.61E+05	1	2	beam-type collision-induced dissociation	60	5952	178.07222	5.41E+04	116.04744	0.82	1, 124
1	6183	7.28	340.150024	7.85E+04	na	2	beam-type collision-induced dissociation	60	6242	224.102753	3.10E+04	116.047272	0.62	1, 124
1	6435	7.57	252.109085	1.04E+05	1	2	beam-type collision-induced dissociation	60	6492	136.061661	3.99E+04	116.047424	0.69	1, 115, 124
1	6794	7.98	228.097778	2.52E+05	na	2	beam-type collision-induced dissociation	60	6847	112.050697	1.78E+04	116.047081	2.27	1, 85
1	7082	8.3	242.113373	1.16E+04	na	2	beam-type collision-induced dissociation	60	7123	126.066032	1.40E+04	116.04734	0.03	1
1	7082	8.31	280.090302	2.50E+05	1	2	beam-type collision-induced dissociation	60	7133	164.043304	7.81E+04	116.046997	2.99	1
1	7082	8.32	290.134613	6.87E+04	na	2	beam-type collision-induced dissociation	60	7143	174.087036	1.72E+04	116.047577	2.01	1, 124, 126
1	7405	8.67	258.108398	8.50E+06	1	2	beam-type collision-induced dissociation	60	7442	142.060913	6.54E+04	116.047485	1.22	1, 124
1	8478	9.93	252.109009	1.86E+05	1	2	beam-type collision-induced dissociation	60	8526	136.061813	4.56E+04	116.047195	1.28	1, 116, 124
1	8478	9.94	229.081787	4.84E+04	na	2	beam-type collision-induced dissociation	60	8537	113.034416	1.74E+04	116.047371	0.23	1, 90
1	8622	10.09	259.092316	6.78E+05	1	2	beam-type collision-induced dissociation	60	8668	143.04483	1.60E+04	116.047485	1.22	1, 85, 101, 124, 125
1	8622	10.1	281.07431	6.80E+05	1	2	beam-type collision-induced dissociation	60	8670	165.027008	3.36E+04	116.047302	0.36	1, 124
1	9269	10.84	284.09906	5.86E+04	na	2	beam-type collision-induced dissociation	60	9310	168.051743	1.33E+04	116.047318	0.23	1, 124, 125
1	9269	10.86	294.06958	6.06E+04	na	2	beam-type collision-induced dissociation	60	9334	178.022522	1.54E+04	116.047058	2.46	1

The output results were automatically classified based on the **tag [1]**.

Ions classified with Tag [1] are accompanied by annotations of all matching tags.

An example of output results from feature annotation and classification (Format-2)

■ The output results can be viewed in excel



When using a tag with **category 3 only** (see [Appendix A4](#))

Information of detected MS ¹ peak pair (MS ¹ -MS ^{1'}) with specific delta mass								Information of detected features			
Classification by tag	MS Level	MS1 Scan Number	RT (min)	Precursor Ion m/z	Precursor Ion Intensity	Precursor Ion' m/z	Precursor Ion' Intensity	Delta Mass (Da)	Mass Error (ppm)	Artificial Adducts	Tags matched
124	1	2269	2.68	268.103972	3.69E+05	152.056664	7.74E+04	116.047308	0.31317	na	2, 83, 124, 125
124	1	2301	2.72	282.119527	2.73E+06	166.072272	3.24E+05	116.047255	0.768416	na	2, 12, 83, 124, 125, 126
124	1	2337	2.76	282.119566	1.61E+07	166.072322	2.26E+05	116.047244	0.862321	M+H	2, 124, 125, 126
124	1	2337	2.76	252.109154	5.18E+04	136.061833	9.45E+05	116.047322	0.193297	M+H	2, 124, 125, 126
124	1	2445	2.88	282.119651	6.97E+06	166.072299	8.55E+04	116.047352	0.072958	na	2, 12, 83, 85, 124, 125, 126
124	1	2481	2.92	282.119589	4.68E+06	166.072272	6.32E+04	116.047316	0.237533	na	2, 124, 125, 126
124	1	2517	2.97	252.109028	3.36E+04	136.061869	6.83E+05	116.047159	1.596199	na	1, 124
124	1	2553	3.01	268.104079	7.91E+04	152.056488	1.43E+04	116.047591	2.126877	na	2, 83, 124, 125
124	1	2589	3.05	282.119682	1.90E+06	166.072339	6.81E+04	116.047343	0.01156	na	2, 12, 83, 124, 125, 126
124	1	2697	3.17	268.104118	7.97E+04	152.056628	1.87E+04	116.04749	1.256413	M+H	2, 83, 124, 125
124	1	2697	3.17	252.108951	3.10E+04	136.061897	5.84E+05	116.047054	2.501978	M+H	2, 83, 124, 125
124	1	2733	3.21	282.11972	8.10E+05	166.072384	4.49E+04	116.047336	0.073181	na	2, 12, 124, 125, 126
124	1	2841	3.34	268.104094	1.05E+05	152.056656	2.03E+04	116.047438	0.812335	na	2, 124, 125
124	1	2877	3.38	282.119548	6.35E+05	166.072311	4.72E+04	116.047237	0.921835	na	2, 12, 124, 125, 126
124	1	2877	3.38	252.108668	3.75E+04	136.06181	5.91E+05	116.046857	4.193028	na	2, 12, 124, 125, 126
124	1	2913	3.42	282.119506	6.63E+05	166.072277	6.56E+04	116.047229	0.987507	na	2, 124, 125, 126
124	1	2985	3.5	268.104049	8.00E+04	152.056569	1.14E+04	116.047481	1.179333	na	2, 124, 125
124	1	3021	3.55	282.119591	3.83E+05	166.072323	4.89E+04	116.047269	0.648275	na	2, 124, 125, 126

An example of output results from feature annotation and classification (Format-3)

■ The output results can be viewed in excel



When using a tag with categories 1 and 3, or 2 and 3 (see [Appendix A4](#))

Information of detected MS ¹ peak pair (MS ¹ -MS ^{1'}) with specific delta mass									Information of detected MS ² ions			Information of detected features				
Classification by tag	MS1 Scan Number	RT (min)	Precursor Ion m/z	Precursor Ion Intensity	Precursor Ion' m/z	Precursor Ion' Intensity	Precursor Charge	Artificial Adducts	Collision Type	Collision Energy	MS2 Scan Number	Product Ion m/z	Product Ion Intensity	Delta Mass (Da)	Mass Error (ppm)	Tags matched
1	2877	3.38	252.108627	172477	na	na	na	na	beam-type collision-induced dissociation	60	2934	136.061523	16350	116.047104	2.06914	1, 124
124	2877	3.38	252.108668	37494	136.06181	590747	na	na	beam-type collision-induced dissociation	60	2934	na	na	116.046857	4.193028	1, 124
1	5897	6.88	294.119659	260874	na	na	1	na	beam-type collision-induced dissociation	60	5952	178.07222	54128	116.04744	0.823588	1, 124
124	5897	6.88	294.119626	200344	178.072346	92776	1	na	beam-type collision-induced dissociation	60	5952	na	na	116.04728	0.5541	1, 124
1	6183	7.21	340.150024	78517	na	na	na	na	beam-type collision-induced dissociation	60	6242	224.102753	31030	116.047272	0.622776	1, 124
124	6183	7.21	340.150061	32461	224.102921	42516	na	na	beam-type collision-induced dissociation	60	6242	na	na	116.04714	1.755252	1, 124
1	6435	7.5	252.109085	103807	na	na	1	na	beam-type collision-induced dissociation	60	6492	136.061661	39865	116.047424	0.6921	1, 115, 124
124	6435	7.5	252.109108	83616	136.061877	2128054	1	na	beam-type collision-induced dissociation	60	6492	na	na	116.047232	0.967764	1, 115, 124
1	7082	8.25	290.134613	68731	na	na	na	na	beam-type collision-induced dissociation	60	7143	174.087036	17248	116.047577	2.006977	1, 124, 126
124	7082	8.25	290.134682	44472	174.087219	13692	na	na	beam-type collision-induced dissociation	60	7143	na	na	116.047463	1.028847	1, 124, 126

3-4. Peak Confirmation via MZmine (optional function) [Step-4; optional]




Fig-12. Interface of “Peak Confirmation via MZmine” in *FeatureHunter*

- 1 Select the “(optional) Peak Confirmation via MZmine”.
- 2 Click “Browse” to import the .tsv file that was generated from Feature Annotation and Classification (Function-2; [step-3](#)).
- 3 Click “Browse” to import the .csv file that was outputted from MZmine¹.
- 4 Click “Add” to submit the data file imported from *FeatureHunter* (.tsv) and MZmine (.csv) for data integration².



- The configuration of the required elements in MZmine's output file is described in [Appendix A6](#).
- *FeatureHunter* enables integrating single, multiple results (from different tags), or aligned results (from different samples) with the MZmine file. (see [Appendix A7](#))

3-4. Peak Confirmation via MZmine (optional function) [Step-4; optional]





Fig-12. Interface of "Peak Confirmation via MZmine" in *FeatureHunter*

- ⑤ **Set the parameters:**
 - (a) Specify a "Project Name" which will be automatically added as a prefix to the output file name. For example: **Integrate**_file name
 - (b) Set the "Mass Tolerance" with ppm for intersecting the precursor ion m/z between *FeatureHunter's* and MZmine's output lists.
 - (c) Set the "RT Tolerance" with seconds for intersecting the precursor ion RT between *FeatureHunter's* and MZmine's output lists.
- ⑥ Specify the output path and folder. (The output files will be in .csv format.)
- ⑦ Click "RUN" to start execution.

Finished



Output success




Fig-13. The integrated results between *FeatureHunter* and MZmine are output in the specified output folder.

An example of output results from data integration

■ The output results can be viewed in excel



Information from
FeatureHunter

Information from
MZmine

Classification by tag	MS1 Scan Number	RT (min)	Precursor Ion m/z	Precursor Ion Intensity	Precursor Ion Charge	MS Level	Collision Type	Collision Energy	MS2 Scan Number	Product Ion m/z	Product Ion Intensity	Delta Mass (Da)	Mass Error (ppm)	Tags matched	MZmine ID	MZmine m/z	MZmine RT (min)	Peak area
1	2255	7.99	494.136383	7.45E+06	1	2	collision-induced dissociation	40	2295	378.088562	1220189	116.047821	4.110778	1, 2, 7, 124,	2938	494.1368	7.905525	6.73E+07
1	2255	8.01	559.147583	2.94E+05	1	2	collision-induced dissociation	40	2308	443.100342	98717	116.047241	0.885751	1, 4	NA	NA	NA	NA
1	2255	8.01	324.129883	2.81E+05	na	2	collision-induced dissociation	40	2311	208.082291	119030	116.047592	2.138464	1, 115, 117	NA	NA	NA	NA
1	2255	8.02	230.11331	2.42E+05	1	2	collision-induced dissociation	40	2323	114.065712	145648	116.047598	2.188721	1	NA	NA	NA	NA
1	2293	8.04	228.097809	9.79E+06	1	2	collision-induced dissociation	40	2336	112.050041	149816	116.047768	3.650572	1, 85, 115	74	228.098	7.946975	1.56E+08
1	2293	8.06	474.182556	9.67E+05	na	2	collision-induced dissociation	40	2362	358.135284	57329	116.047272	0.622776	1, 124	NA	NA	NA	NA
1	2331	8.11	562.111755	1.69E+05	1	2	collision-induced dissociation	40	2400	446.063873	51605	116.047882	4.636729	1	3682	562.111	7.905525	2.65E+06
1	2369	8.11	344.145508	1.41E+05	na	2	collision-induced dissociation	40	2405	228.09761	22417	116.047897	4.768217	1, 4, 85, 124	NA	NA	NA	NA
1	2369	8.12	455.187836	1.94E+08	1	2	collision-induced dissociation	40	2406	339.140961	42771	116.046875	4.041454	1, 12	2678	455.1884	7.905525	3.07E+09
1	2369	8.15	554.080566	2.14E+05	na	2	collision-induced dissociation	40	2434	438.033112	73726	116.047455	0.955076	1	NA	NA	NA	NA
1	2404	8.17	288.041168	7.21E+05	na	2	collision-induced dissociation	40	2451	171.994156	10244	116.047012	2.858065	1	883	288.0413	7.988232	1.06E+07
1	2404	8.19	282.149597	9.25E+04	na	2	collision-induced dissociation	40	2474	166.102386	13120	116.047211	1.148726	1	NA	NA	NA	NA
1	2441	8.2	310.088409	9.33E+06	1	2	collision-induced dissociation	40	2483	194.040573	818125	116.047836	4.242266	1	1167	310.0885	7.988232	1.13E+08
1	2477	8.28	296.072845	4.07E+05	na	2	collision-induced dissociation	40	2549	180.02507	25570	116.047775	3.716316	1	964	296.0727	7.988232	4.04E+06
1	2550	8.34	290.176208	1.47E+05	1	2	collision-induced dissociation	40	2604	174.128403	26686	116.047806	3.979291	1	NA	NA	NA	NA
1	2588	8.36	261.127411	1.41E+06	na	2	collision-induced dissociation	40	2627	145.079529	1520372	116.047882	4.636729	1	512	261.1274	8.277544	1.48E+07
1	2626	8.41	343.117889	5.03E+05	1	2	collision-induced dissociation	40	2673	227.070267	37964	116.047623	2.401439	1, 124	1572	343.1178	8.277544	3.89E+06
1	2626	8.42	258.108643	4.96E+04	na	2	collision-induced dissociation	40	2683	142.06073	59294	116.047913	4.899704	1, 2	436	258.1086	9.531694	7.58E+07



NA indicates that the ions were not detected in MZmine.
It is recommended to exclude these ions from the final results due to their low peak quality and confidence.

4. Appendix-A1

Visual Studio 2015, 2017, 2019, and 2022

This table lists the latest supported English (en-US) Microsoft Visual C++ Redistributable packages for Visual Studio 2015, 2017, 2019, and 2022. The latest supported version has the most recently implemented C++ features, security, reliability, and performance improvements. It also includes the latest C++ standard language and library standards conformance updates. We recommend that you install this version for all applications created using Visual Studio 2015, 2017, 2019, or 2022.

Unlike older versions of Visual Studio, which have infrequent redistributable updates, the version number isn't listed in the following table for Visual Studio 2015-2022 because the redistributable is updated frequently. To find the version number of the latest redistributable, download the redistributable you're interested in using one of the following links. Then, look at its properties using Windows File Explorer. In the **Details** pane, the **File version** contains the version of the redistributable.

Architecture	Link	Notes
ARM64	https://aka.ms/vs/17/release/vc_redist.arm64.exe	Permalink for latest supported ARM64 version
X86	https://aka.ms/vs/17/release/vc_redist.x86.exe	Permalink for latest supported x86 version
X64	https://aka.ms/vs/17/release/vc_redist.x64.exe	Permalink for latest supported x64 version. The X64 Redistributable package contains both ARM64 and X64 binaries. This package makes it easy to install required Visual C++ ARM64 binaries when the X64 Redistributable is installed on an ARM64 device.

Download and install

Fig-A1. The page for downloading Visual C++ redistributable package

Visual C++ redistributable package can be downloaded from the official website of Microsoft:
<https://learn.microsoft.com/en-us/cpp/windows/latest-supported-vc-redist?view=msvc-170>

4. Appendix-A2





Fig-A2. The page for downloading .NET Framework

.NET Framwork can be downloaded from the official website of Microsoft:

<https://dotnet.microsoft.com/en-us/download/dotnet-framework/net472>

4. Appendix-A3



MSConvert can be downloaded from the website:

<https://proteowizard.sourceforge.io/download.html>

Fig-A3. Recommended setting parameter for converting data from vendor's format to mzML

4. Appendix-A4

	A 1	B 2	C 3	D 4	E 5	F 6	G 7	H 8	I 9
	Tag No.	Tag Category	Tag Parameters (Da or m/z)	Mass Tolerance (ppm)	Intensity Cutoff (height)	Top N Ion by Intensity (MS2; only for tag category 1-2)	Consecutive Data Points (MS1; only for tag category 3)	Minimum Intensity Ratio (MS1/MS1'; only for tag category 3)	Maximum Intensity Ratio (MS1/MS1'; only for tag category 3)
1									
2	1	1	116.047344	5	10000	20	0	0	0
3	2	1	132.042258	5	10000	20	0	0	0
4	3	1	146.057908	5	10000	20	0	0	0
79	78	2	151.049409	5	10000	20	0	0	0
80	79	2	152.056686	5	10000	20	0	0	0
81	80	2	135.030137	5	10000	20	0	0	0
106	105	3	3.01883	5	10000	0	5	0.7	1.3
107	106	3	-3.01883	5	10000	0	5	0.7	1.3

Fig-A4. The management of Feature List.csv

1. Tag No.: Tag No. is utilized to characterize the tag parameters in *FeatureHunter*. Each tag represents a specific feature or characteristic that can be assigned to annotated data for classification purpose. The default Tag No. ranges from 1 to 132. Users have the flexibility to manage the feature list by adding, editing, or removing tags as needed.


2. Tag Category.: Tag category is utilized to classify tags into three types.

These categories include:

Tag category 1: neutral loss of a particular functional group.

Tag category 2: specific diagnostic product ion m/z in MS^2 .

Tag category 3: a fixed delta mass in MS^1 .



3. Tag Parameters: Tag parameters are utilized to extract the diagnostic features of the analytes of interest obtained from MS^1 or MS^2 raw data.

4. Mass Tolerance: Mass tolerance is utilized for tag parameters with a default value of 5 ppm. Users can adjust this value to meet their data quality requirements.

5. Intensity Cutoff (height): Intensity cutoff is utilized to filter features (such as product ions or peak pairs) based on their peak height, retaining only those above the cutoff value. The default value is 10000. Users can adjust this value to meet their data quality requirements.

6. Top N ion by intensity: The Top N ions filter is utilized to retain product ions with intensity rankings within the top N in the MS² spectrum. This function is applicable only to categories 1 and 2. The default value is 20. Users can adjust this value to meet their data quality requirements. The function can be disabled by setting the value to zero.

7. Consecutive Data Points: Consecutive data points are utilized to filter specific MS¹ peak-pairs with a fixed delta mass across consecutive MS¹ spectra (data points). This function is applicable only to category 3. The default value is 5. Users can adjust this value to meet their data quality requirements. The function can be disabled by setting the value to zero.

8. Minimum Intensity Ratio: The minimum intensity ratio is utilized to filter specific MS¹ peak-pairs based on a specified intensity ratio. The default value is 0.7. Users can adjust this value to meet their data quality requirements. The function can be disabled by setting the value to zero.

- Intensity ratio of peak-pairs:
$$\frac{\text{Precursor ion}}{\text{Precursor ion}'}$$

9. Maximum Intensity Ratio : The maximum intensity ratio is utilized to filter specific MS¹ peak-pairs based on a specified intensity ratio. The default value is 1.3. Users can adjust this value to meet their data quality requirements. The function can be disabled by setting the value to zero.

- Intensity ratio of peak-pairs:
$$\frac{\text{Precursor ion}}{\text{Precursor ion}'}$$

4. Appendix-A5

	①	②	③
	A	B	C
	Artificial Adduct No.	Artificial Adduct m/z	Artificial Adduct Name
1			
2	1	21.981945	M+Na-H
3	2	37.955882	M+K-H
4	3	17.026549	M+NH3
5	4	32.026214	M+MeOH
6	5	60.057514	M+IPA
7	6	78.013935	M+DMSO
8	7	46.041864	M+EtOH


Fig-A5. The management of Artificial Adducts List

1. Artificial Adduct No.: The ID of each artificial adduct in the list.


2. Artificial Adduct m/z: Artificial Adduct m/z is used to annotate the adduct ion of the precursor ion in the MS¹ spectrum. The default artificial adduct list comprises 13 common artificial adducts. Users have the flexibility to manage the feature list by adding, editing, or removing artificial adducts as needed.

3. Artificial Adduct Name: Artificial Adduct name is used for characterizing the identity of artificial adducts.


5. Appendix-A6



Intersection
 $[A \cap B]$



Union
 $[A \cup B]$



Difference
 $[A - B]$

Example:

■ $[1]:$

Intersection: 1 Add

■ $[1 \cap 124]:$

Intersection: 1,124 Add

■ $[2 \cup 3 \cup 4]:$

Union: 2,3,4 Add

■ $[1 \cap 124] - [2]:$

Intersection: 1,124
Difference: 2 Add

■ $[2 \cup 3] - [4]:$

Union: 2,3
Difference: 4 Add

■ $[1] \cap [2 \cup 3] \cap [7 \cup 8]:$


Step-1 Intersection: 1
Union: 2,3 Add

Step-2 Union: 7,8 Add

■ $[42] \cap [58] \cap [78 \cup 79 \cup 80] - [115 \cup 116]:$

Intersection: 42,58
Union: 78,79,80
Difference: 115,116 Add

5. Appendix-A7



The screenshot shows the 'Please set the parameters' dialog box in FeatureHunter. The window is titled 'Please set the parameters' and has a standard macOS-style title bar with red, yellow, and green buttons. The interface is organized into several sections on the left: 'Feature lists', 'Filename', 'Field separator', 'Export common elements', 'Export data file elements', 'Export quantitation results and other information', 'Identification separator', and 'Filter rows'. The 'Feature lists' section shows '0 selected' and a dropdown menu set to 'As selected in main window'. The 'Filename' section shows the file path '/20231115_QC_aligned.csv' with 'last' and '...' buttons. The 'Field separator' section has a text input field containing a comma ','. The 'Export common elements' section contains a list of checkboxes: 'Export row ID' (checked), 'Export row m/z' (checked), 'Export row retention time' (checked), 'Export row identity (main ID)' (unchecked), 'Export row identity (all IDs)' (unchecked), 'Export row identity (main ID + details)' (unchecked), and 'Export row comment' (unchecked). To the right of this list are 'All' and 'Clear' buttons. The 'Export data file elements' section contains another list of checkboxes: 'Peak RT start' (unchecked), 'Peak RT end' (unchecked), 'Peak duration time' (unchecked), 'Peak height' (unchecked), 'Peak area' (checked), 'Peak charge' (unchecked), and 'Peak # data points' (unchecked). To the right of this list are 'All' and 'Clear' buttons. The 'Export quantitation results and other information' section has an unchecked checkbox. The 'Identification separator' section has a text input field containing a semicolon ';'. The 'Filter rows' section has a dropdown menu set to 'ALL'. At the bottom of the window are 'OK', 'Cancel', and 'Help' buttons.

FeatureHunter has been verified to support data integration with file exported from both older (2.53) and latest versions (3.90) of MZmine.

The processed results, exported as a .CSV file from MZmine, should include the following elements:



- (1) Row ID
- (2) Row m/z
- (3) Row retention time
- (4) Peak area



When exporting data from MZmine 3, select "Export CSV (legacy MZmine 2)".

Fig-A6. The configuration of required elements for data output from MZmine

5. Appendix-A8



Report Issues



We welcome your feedback and comments, as they can greatly contribute to enhancing the development of the software.

If you encounter any issues with *FeatureHunter*, please include the following information in your issue report:

- (1) *The version of FeatureHunter.*
- (2) *Screenshots of error messages or unexpected results.*
- (3) *Any additional information that could help us better understand the challenges or bugs you're encountering.*

Please leave your feedback and comments at the following link:

https://msomics.abrc.sinica.edu.tw/FeatureHunter/?page_id=271



Contact Us



Agricultural Biotechnology Research Center,
Academia Sinica, Taiwan



Exposomics Research Laboratory,
Chung Shan Medical University, Taiwan

Dr. Yet-Ran Chen

E-mail: yetran@gate.sinica.edu.tw

Mr. Wei-Hung Chang

E-mail: whchang@gate.sinica.edu.tw

Dr. Mu-Rong Chao

E-mail: mrchao@csmu.edu.tw

Dr. Chiung-Wen Hu

E-mail: windyhu@csmu.edu.tw

Dr. Yuan-Jhe Chang

E-mail: jeffchang@csmu.edu.tw